

Double-spot LED Matrix Light



with electronic dimmer

Quality of Inspection

Double-spot LED illuminator containing 2x18 ultra-bright white LED lamps mounted on gooseneck arm for various oblique incident light applications. Equipped with a universal clamp compatible with all Inspectis microscope Stand models.



For high-contrast illumination of different specimens and applications.

Spot-light is powered directly via AUX output of Inspectis Digital Microscopes for excellent cable management. ON-OFF switch and brightness control is integrated. Extra power supply for compatibility to other devices are provided as option.

Specification:

Number of light elements	2x18 pcs \varnothing 5 mm LED lamps
Colour temperature	8600 K
Size of LED panel	78x32.5x23 mm
Power	12VDC, 600mA via AUX or External power supply
Length of gooseneck arm	280mm
Total weight	480 g. max

Ordering Number/Name

HD-126 Inspectis Double-Spot LED Matrix Light
HD-024-S External Power Supply (optional)



Reseller:

V1, 2018

www.inspect-is.com